

CASE REPORT



Treatment of Focal Arterial Embolus With Stand-Alone Mechanical Thrombectomy



John A. Phillips, MD
OhioHealth Vascular
Institute/Riverside
Methodist Hospital
Columbus, OH

PATIENT PRESENTATION

An 82-year-old woman presented with symptoms of acute-onset pain and paresthesia. The patient's initial vascular exam was abnormal. The patient had a complex prior medical history, including chronic kidney disease, atrial fibrillation, heart failure with preserved ejection fraction, type 2 diabetes, hypertension, and dyslipidemia.

DIAGNOSTIC FINDINGS

After the abnormal vascular exam, an angiogram revealed a focal embolic occlusion of the right common femoral artery (CFA) (**Figure One**). A noninvasive exam showed an ankle-brachial index (ABI) of 0 on the right leg and 1.26 on the left leg. A 5 Fr, 10 cm Pinnacle® Destination® Sheath was placed, and an arteriogram was obtained, which confirmed an atrial fibrillation embolus in the right leg CFA.

Pre-procedure Angiogram
of the Right CFA



Figure One

The Pounce™ System Basket
Retracted into Funnel in the
Right CFA




Figure Two

TREATMENT

The patient was immediately started on heparin intravenously to prepare for an intervention. Left groin access was obtained by micropuncture, and a 7 Fr, 45 cm Flexor® Ansel Guiding Sheath was placed contralaterally into the right CFA.

The Pounce™ Thrombectomy System was prepared and a .035 Magic Torque™ Guidewire was used to cross the lesion. The Pounce system delivery catheter was placed distal to the embolus, the Pounce system basket wire was delivered through the delivery catheter into the mid superficial femoral artery (SFA), and the Pounce system funnel catheter was advanced and deployed within the CFA with some draping into the bifurcation profunda SFA to clear the distal CFA.

The baskets were then pulled back into the funnel (**Figure Two**), capturing the embolus, and the device was removed through the guiding sheath. After cleaning the Pounce system, the baskets and funnel were deployed again in similar

CASE REPORT

positions and another pass was made. Another angiogram was obtained (**Figure Three**), showing complete removal of the embolus in the two passes (**Figure Four**). No ancillary treatments (additional thrombectomy device, drug-coated balloon, or stent) were used.



Figure Three




Figure Four

POST PROCEDURE OUTCOME

The patient was transferred back to the hospital floor symptom-free. Noninvasive studies performed after the procedure showed improvement of the ABI to 1.05 on the right leg and 1.36 on the left leg.

The Pounce™ Thrombectomy System provided prompt, on-table restoration of arterial flow for a patient with an embolic occlusion of the right CFA. No thrombolytics or other adjunctive therapies were used during the procedure.

I always worry about treating a thrombus in patients with underlying peripheral artery disease because their arteries may not be healthy. Having said that, I found the Pounce system to be fairly atraumatic. With this device, I believe we now have an endovascular option for arterial thromboembolic phenomena that previously required surgical intervention.



Surmodics, Inc.
7905 Golden Triangle Drive, Suite 190 Eden Prairie, MN 55344 USA
Toll Free: 888-626-8501 Phone: 952-500-7400
www.pouncesystem.com

Caution: Federal (US) law restricts this device to sale by or on the order of a physician. Please refer to Instructions for Use for indications, contraindications, warnings, and precautions.

SURMODICS, POUNCE, and SURMODICS and POUNCE logos are trademarks of Surmodics, Inc. and/or its affiliates. Third-party trademarks are the property of their respective owners. All images used with permission of author.

©2023 Surmodics, Inc. All rights reserved. SRDX-PTS-764 A 02/23